How to Archive a Project?

Droplet Lab Projects | Surface Science Instruments & Solutions · https://dropletlab.com/knowledge-doc/how-to-archive-a-project/

1. From the Projects screen, long press a project you want to archive. Doing so selects the project and a checkbox appears at the right side displaying that the project has been selected.



2. Upon selection, tap the Archive Project button displayed at the top-right corner of the screen.


3:55
31%
Projects
Archive Projects
TEST PROJECT
Created: 16 Jul 2023 3:55 am
Add new project