Ability to Back Up App Data

Release Notes · https://dropletlab.com/knowledge-doc/ability-to-back-up-app-data/

With this enhancement, the user has the ability to back up all application data. This adds to the convenience of the user in case they lose access to their account.

How to Backup App Data?

1. Upon successful login, click the hamburger icon.
2. From the screen displayed, click Back Up. This starts a process to save the data.
3. Once the process completes, a notification shows the destination of the file.
4. Connect the phone to a computer to copy the data or use any cloud storage app to upload it to your cloud.