

App does not Start. What Should I Do?

Troubleshooting · https://dropletlab.com/knowledge-doc/app-does-not-start-what-should-i-do/


1. In case the app does not start, tap the Settings icon from your phone screen.




2. From the Settings app, tap Apps. This opens the Apps screen.




3. Here, locate and tap the Droplet Lab app. Doing so takes you to the App details screen.



4. Tap the Permissions button.



5. From the App Permissions screen, make sure that all permissions are allowed.


11:43
61%
App permissions
Energy
Allowed
Camera
Contacts
Files
Location
Music and audio
Nearby devices
Notifications
Phone
Photos and videos