

How to Check for App Version and Software Updates?

App Installation · https://dropletlab.com/knowledge-doc/how-to-check-for-app-version-and-software-updates/


1. From the app’s landing screen, tap the three-line icon located at the bottom-left corner of the screen.



2. From the options displayed, click About. Doing so reveals a popup displaying the current version of the application.



3. To check for updates, click the Check for Updates button.


6:34
Settings
About
Tilted
Version: 1.4.0.10
ACCEPT
CHECK FOR UPDATES
Report an issue
Show help Messages